

Strzykawki do buforów, modyfikowane, uszczelka PTFE



Strzykawki do buforów, modyfikowane, uszczelka PTFE

Charakterystyka

Strzykawki do buforów, modyfikowane, uszczelka PTFE

W zapytaniu ofertowym należy podać parametry podane w tabeli

Volume	AD	DAD	DX	C	CX	TLL	TLL-X	RQ	XL
25 µl	-	-	2627517	2606517	2606527	2641517	2643517	2644517	2645517
50 µl	2605617	2626617	2627617	2606617	2606627	2641617	2643617	2644617	2645617
100 µl	2605717	2626717	2627717	2606717	2606727	2641717	2643717	2644717	2645717
250 µl	2605817	2626817	2627817	2606817	2606827	2641817	2643817	2644817	2645817
500 µl	2605917	2626917	2627917	2606917	2606927	2641917	2643917	2644917	2645917
1.0 ml	2605017	2626017	2627017	2606017	-	2641017	2643017	2644017	2645017
2.5 ml	2605037	2626037	-	2606037	-	2641037	2643037	2644037	2645037
5.0 ml	2605057	-	-	2606057	-	2641057	2643057	2644057	2645057
10.0 ml	2605077	-	-	2606077	-	2641077	2643077	2644077	2645077
Front Fitting	M8x0,75	M8x0,75	Luer	1/4"-28 UNF	1/4"-28 UNF	Luer	Luer	Luer	Luer
Side Connection	-	M6x1	M6x1	-	-	-	-	-	-
Plunger Top	-	6-32 UNC	6-32 UNC	-	6-32 UNC	-	6-32 UNC	+ Screw	-

Volume	H-PD	XE	XP
12.5 µl	2666415	-	-
25 µl	2666515	-	2647517
50 µl	2666615	-	2647617
100 µl	2666715	-	2647717
125 µl	2666714	-	-
250 µl	2666815	-	2647817
500 µl	2666915	2646927	2647917
1.0 ml	2666014	2646027	2647017
1.25 ml	2666015	-	-

2.5 ml	2666035	2646047	2647037
5.0 ml	2666055	2646067	2647077
12.5 ml	2666074	-	-
Front Fitting	Luer	1/4"x28 UNF	1/4"x28 UNF
Plunger Top	6-32 UNC	-	-

Dostępne modele

Model

Opis

Do buforów

Do wyboru parametry podane w tabeli